

Art. 8706 Bain Marie system stainless steel

Date last modified: 5-8-2022, Changes reserved for printing- and typing errors.

Product image



Product information

The standard melting tank works according to the Au-Bain-Marie system. This means that the product that has to be heated is heated in a can or tray with hot water. The water in this tank is heated by a thermostat-controlled element maintained at a variably set temperature of 30-85°C.

Installation

1. Always place the tank on a stable surface.
2. Make sure the faucet is in the closed position.
3. Place the stainless steel cans/trays.
4. Connect the device to an earthed socket with 230 Volt.
5. Check whether the device is still switched off. The thermostat must therefore be set to zero.
6. The tank is ready for use

Use of appliance

1. Fill the tank with water until the filled cans and/or trays (paraffin) float.
2. Check if there is enough water in the tank. The more water there is in the tank, the better the heat transfer to melt the paraffin.
3. Switch on the thermostat to the desired temperature (approx. 70-80°C) and wait until the paraffin has completely melted.
4. The device is ready for use.

Note!

Regularly check the amount of water in the tank and prevent it from boiling dry! The more water, the better the heat transfer. For long intervals of use, set the temperature back to 60°C. This is to maintain the quality of the paraffin and to save the tank. Always switch off the device after use. Also remove the plug from the socket. This prevents anything from going wrong at all times. Never use the melting tank for improper use. This is to prevent damage and increase the lifespan. Improper use also invalidates the warranty.

Cleaning

1. Switch off the appliance and remove the plug from the socket.
2. To clean the cans/trays, it is wise to have them as empty as possible.
3. Pour the last remaining paraffin from the cans/trays and store it in a disposable container that is heat-resistant. This (solidified) remainder can be used again later on in the cans.
4. The cans and trays are easy to clean with cleaning paper when they are warm.
5. When the cans/trays are clean, they can be cleaned extra with some soapy water.
6. Rinse the soapy water well and dry everything thoroughly before using it again.
7. The tank can be emptied and scrubbed after removing the dipping cans.
8. The tank is then ready for use again.

Specifications

Length	540mm
Width	340mm
Height	240mm
Depth	165mm
Wattage	1200Watt
Voltage	230Volt
Thermostat setting	30-85°C
Variable layout	✓
Maximum security	✓

Warranty


The melting tank comes with a one-year warranty against material and construction defects, provided that the device has been used and maintained in accordance with the instructions. Also, the melting chain must not be mishandled or abused in any way. This only applies to parts, not to labor and transport. In the event of malfunctions, please contact your supplier in order to have the malfunction resolved with the original parts.

Warning: as this involves heating a liquid, operation of the melting tank should only be done by qualified personnel.

The device complies with EU directives 73/23 EEC and 89/336 EEC, all technical changes reserved.



Electric scheme



Always read the manual before using the device